

# RT300 Supine

## Early mobility cycling



# RT300 Supine

## Early mobility cycling

Unique design makes RT300 Supine easy to use. RT300 is mounted on an expandable table fitting over and around beds allowing RT300 to be brought to the patient. RT300's height is easily adjusted up to 46.5" (118cm) via the motorized table top.

- Accommodates beds up to 46" (116 cm) wide including bariatric beds
- Anti-microbial touch screen and cleanable surfaces
- Adult and pediatric pedals available
- Adjustable cranks cater to patient size and asymmetry
- Cleanable gel pads for skin protection

### Internet connectivity with HIPAA compliant database



Facilitates automatic download of patient settings and uploading of performance data making RT300 Supine attractive for clinical application and research.

### Arm cycle

Switch out pedals for hand grips and arm cycle while supine, seated, or standing.

### Available with optional integrated functional electrical stimulation (FES)

When using the integrated FES, the following benefits can be achieved even for patients who are very weak or cannot participate volitionally.

- Prevent or retard disuse atrophy
- Increase local blood circulation
- Maintain or increase range of motion
- Relax muscle spasms



**Enhance your early mobility program  
with RT300 Supine.**

For more information, contact us at  
**supine@RT300.com**  
or call **1-800-609-9166, x302.**

